

ADVENTURER SEMI-MICRO AX85 – 30852308

SPECIFICATIONS

Maximum Capacity

82 g

Readability

0.00001 g

Pan Size

80 mm

Internal Calibration

AutoCal™ - Automatic

Draftshield

Included

Dimensions

9.05 in (354 mm x 340 mm x 230 mm)

Legal for Trade

Not Applicable

Auxiliary Display Model

Available as an Accessory

Auto Door	Not Applicable
Working Environment	10°C – 30°C, 80%RH, non-condensing
Net Weight	5.1 kg
In-Use Cover	Included
Communication	USB; RS232 (Included); USB Host
Power	AC Adapter (Included)
Units of Measurement	Singapore Tael; Ounce Troy; Pennyweight; Grain; Tical; Custom; Milligram; Momme; Newton; Baht; Taiwan Tael; Gram; Hong Kong Tael; Pound; Tola; Mesghal; Carat; Ounce
Display	Full-Color 4.3 in Touchscreen WQVGA
Pan Construction	Stainless Steel

Built-in Ionizer	No
Stabilization Time	8 s
Repeatability, typical	0.00002 g
Minimum Weight (USP, 0.1%, typical)	20 mg
Linearity \pm	0.1 mg
