

COMPASS CX PORTABLE BALANCES – CX1201 – 30428210

SPECIFICATIONS

Maximum Capacity	1,200 g
Readability	0.1 g
Pan Size	5.6 in x 5 in (142 mm x 128 mm)
Battery Life	1001 Hours with Disposable Batteries
Units of Measurement	Gram; Pound:Ounce; Ounce; Newton
Legal for Trade	Not Applicable
Dimensions	1.6 in x 7.9 in x 5.3 in (41 mm x 201 mm x 135 mm) (HxLxW)
Display	LCD Display with white backlight

Working Environment	64°F – 77°F, 80%RH, non-condensing (18°C – 25°C, 80%RH, non-condensing)
Tare Range	To capacity by subtraction
Pan Construction	Stainless Steel
Stabilization Time	2 s
Net Weight	1.1 lb (0.5 kg)
In-Use Cover	No
Transportation Case	Available as an Accessory
Power	3 AA (LR6) Batteries (Included); AC Adapter (Not Included)
Communication	Not Applicable
Auxiliary Display Model	Not Applicable
