

ENDEAVOUR 5000 SERIES E-E51LD0403 EU 30770915

SPECIFICATIONS

Motion	Orbital, 3 mm
Capacity	4.0 kg, up to 3 x 1 L flasks
Speed Range	100 rpm – 800 rpm
Tray Size	11.75 in x 8.75 in (300 mm x 222 mm) (LxW)
Timer	1 minute – 99 hours, 59 minutes
Control	Digital
Speed Accuracy	±2 rpm, Above 100 rpm ±2% of set speed
Drive System	Triple Eccentric; Brushless DC Motor

Tray Construction	Aluminum
Dimensions	16.5 in x 4.7 in x 10.3 in (420 mm x 119 mm x 262 mm) (LxHxW)
Orbit	0.1 in (3 mm)
Display	Custom Back-lit LCD
Working Environment	-10°C – 40°C, 80%RH, non-condensing (14°F – 104°F, 80%RH, non-condensing)
Speed Range	100 rpm - 800 rpm
Net Weight	16.5 lb (7.5 kg)
Safety Certification	CE; TUV
Power	24V DC, 1.25A
Communication	RS232

Width (Shipping)	40 cm
Height (Shipping)	32 cm
Length (Shipping)	56 cm
Gross Weight (Shipping)	10.2 kg
