

# CENTRIFUGE FRONTIER 5000 – FC5306 – 30134156

## SPECIFICATIONS

Speed Range	0 rpm – 6,000 rpm
Maximum Relative Centrifuge Force (x g)	2,000 g
Maximum Capacity(Rotor)	8 x 1.5 / 2.0 ml
Program memories	0
Rotor In Package	Angle rotor 4 x 8 place PCR Strips; Angle 8 x 1.5 ml / 2.0 ml
Acceleration Gear	1
Deceleration Gear	2
Dimensions	6.6 in x 5 in x 5.9 in (168 mm x 128 mm x 150 mm) (LxHxW)

---

Display	No display
Lid lock	Mechanical
Net Weight	2.3 lb (1.1 kg)
Noise Level	≤60±2 dB(A)
Power	100 - 240 V
Pre-cool time	Not applicable
Rotor Compatibility	Angle 8 x 1.5 ml / 2.0 ml; Angle 4 x 8 place PCR Strips
Rotor Recognition Method	Not Applicable
Running Time	Continuous
Safety Certification	CE; TUV; FCC

---

---

Working Environment

41°F – 104°F, 80% RH, Non-condensing (5°C – 40°C, 80% RH, Non-condensing)

---

