

CENTRIFUGE 5000 SERIES MULTI PRO 100-230V – FC5707 – R09 83041073 SPECIFICATION

Speed Range	200 rpm – 6,800 rpm
Maximum Relative Centrifuge Force (x g)	4,445 g
Maximum Capacity(Rotor)	4x50ml RB/FA or 4x15ml RB/FA
Program memories	99
Rotor In Package	83041009 - Angle rotor 4x50ml RB/FA or 4x15ml RB/FA
Acceleration Gear	10
Deceleration Gear	10
Dimensions	14 in x 9.4 in x 10.9 in (356 mm x 236 mm x 277 mm) (LxHxW)

Display	LCD with backlight
Lid lock	Mechanical
Net Weight	22.1 lb (10 kg)
Noise Level	≤60±2 dB(A)
Power	AC Power (Included)
Power Consumption	50 W
Pre-cool time	Not applicable
Rotor Compatibility	Angle rotor 4x50ml RB/FA or 4x15ml RB/FA
Rotor Recognition Method	Manual
Running Time	10 sec to 99 h 59 min 59 sec

Safety Certification

CE; FCC

Working Environment

36°F – 95°F, 80%RH, non-condensing (2°C – 35°C, 80%RH, non-condensing)

