

FRONTIER SERIES 5000 MULTI – FC5706 – 30332131

SPECIFICATIONS

Speed Range	200 rpm – 6,000 rpm
Maximum Relative Centrifuge Force (x g)	4,427 g
Maximum Capacity(Rotor)	6 x 50 ml
Rotor In Package	30130877 - Angle rotor 12 x 15 ml
Acceleration Gear	10
Deceleration Gear	10
Dimensions	15.4 in x 10.5 in x 11.5 in (392 mm x 266 mm x 291 mm) (LxHxW)
Display	LCD with backlight

Lid lock	Motor driven
Net Weight	23 lb (10.5 kg)
Noise Level	≤60±2 dB(A)
Power	AC Power (Included)
Power Consumption	100 W
Pre-cool time	Not applicable
Rotor Compatibility	Swing Out Rotor 6 x 5 ml; Angle 12 x 15ml; Angle 6 x 50 ml
Rotor Recognition Method	Manual
Running Time	Continuous; 10 sec to 99 h 59 min 59 sec
Safety Certification	CE; FCC

Working Environment

36°F – 95°F, 80%RH, non-condensing (2°C – 35°C, 80%RH, non-condensing)

