

CENTRIFUGE 5000 SERIES MULTI PRO 230V – FC5714 – 30314810

SPECIFICATION

Speed Range	200 rpm – 14,000 rpm
Maximum Relative Centrifuge Force (x g)	18,624 g
Maximum Capacity(Rotor)	4 x 100 ml
Program memories	99
Rotor In Package	Not applicable
Acceleration Gear	10
Deceleration Gear	10
Dimensions	19.41 in x 12.99 in x 14.25 in (493 mm x 330 mm x 362 mm) (LxHxW)

Display	LCD with backlight
Lid lock	Mechanical
Net Weight	66 lb (30 kg)
Noise Level	≤63±2 dB(A)
Power	AC Power (Included)
Power Consumption	240 W
Pre-cool time	Not applicable
Rotor Compatibility	Angle 12 x 15 ml RB/FA ID; Angle 24 x 1.5 / 2.0 ml ID Sealable; Angle 30 x 15 ml RB/FA ID; Swing out 4 x 100 ml ID; Angle 30 x 15 ml RB/FA ID; Angle 24 x 1.5 / 2.0 ml ID Sealable; Swing out 4 x 100 ml ID Sealable; Angle 24 x 1.5 / 2.0 ml ID; Angle 6 x 50 ml RB/FA ID; Swing out 2 x 30 ml ID Sealable
Rotor Recognition Method	Auto

Running Time	Continuous; Quickspin; 10 sec to 99 h 59 min 59 sec
Safety Certification	CE; FCC
Working Environment	36°F – 95°F, 80%RH, non-condensing (2°C – 35°C, 80%RH, non-condensing)

