

CENTRIFUGE 5000 SERIES MULTI PRO 230V – FC5816 – 30314816

SPECIFICATION

Speed Range	200 rpm – 15,000 rpm
Maximum Relative Centrifuge Force (x g)	21,379 g
Maximum Capacity(Rotor)	6 x 250 ml
Program memories	99
Rotor In Package	Not applicable
Acceleration Gear	10
Deceleration Gear	10

Dimensions	21.18 in x 13.94 in x 17.56 in (538 mm x 354 mm x 446 mm) (LxHxW)
Display	LCD with backlight
Lid lock	Motor driven
Net Weight	115 lb (52 kg)
Noise Level	≤61±2 dB(A)
Power	AC Power (Included)
Power Consumption	530 W
Pre-cool time	Not applicable
Rotor Compatibility	Angle 4 x 85ml RB ID Hi; Angle 6 x 250 ml FB ID; Angle 44 x 1.5 / 2.0 ml ID V2; Angle 30 x 15 ml RB/FA ID; Angle 6 x 50 ml RB ID Hi; Rotor Angle 24x1.5/2.0ml ID BIOSEALS V1; Angle 10 x 50 ml FA ID; Angle 20 x 10 ml RB ID Hi; Angle 6 x 85 ml RB ID; Swing out 4 x 250 ml ID; Swing out 2 x 3MTP ID

Rotor Recognition Method	Auto
Running Time	Continous; Quickspin; 10 sec to 99 h 59 min 59 sec
Safety Certification	CE; FCC
Working Environment	36°F – 95°F, 80%RH, non-condensing (2°C – 35°C, 80%RH, non-condensing)
