

GUARDIAN-5000-E-G51HS07C-HOTPLATE & STIRRERS -30500523

SPECIFICATION

Function	Heating and Stirring
Capacity	15 L
Plate Construction	Ceramic
Plate Size	7 in x 7 in (178 mm x 178 mm)
Temp. Range	Ambient +5°C – 500°C
Speed Range	60 rpm – 1,600 rpm
Calibration	Temperature Calibration 3 Plate/3 Probe
Control	Digital

Dimensions	12.2 in x 4.8 in x 8.8 in (309 mm x 122 mm x 223 mm) (LxHxW)
Display	Custom Back-lit LCD
In-Use Cover	Available as an Accessory
Net Weight	5.7 lb (2.6 kg)
Power	230V, 6.0A, 50/60Hz
Power Consumption	1,242 W
Safety Certification	TUV
Speed Range	60 to 1600 rpm
Speed Stability	+/-2%
Temperature Probe	Available as an Accessory

Temp. Stability

+/-3% Plate, below 100°C +/-2°C and +/-1% Probe, below 100°C +/-1°C

Working Environment

41°F – 104°F, 80% RH, Non-condensing (5°C – 40°C, 80% RH, Non-condensing)

