

GUARDIAN 3000 HOTPLATES AND STIRRERS – E-G31HS04C 30680280

SPECIFICATION

Function	Heating and Stirring
Capacity	15 L
Plate Construction	Ceramic
Plate Size	4 in x 4 in (102 mm x 102 mm)
Temp. Range	Ambient +5°C – 500°C
Control	Digital
Dimensions	10.4 in x 4.3 in x 6.4 in (264 mm x 110 mm x 163 mm) (LxHxW)

Power Consumption	550 W
Display	Custom Back-lit LCD
Working Environment	5°C – 40°C, 80% RH, Non-condensing (41°F – 104°F, 80% RH, Non-condensing)
Speed Range	80 to 1600 rpm
Net Weight	5.5 lb (2.5 kg)
Safety Certification	TUV
Temp. Stability	+/-3% Plate, Below 100°C +/-2°C and +/-2% Probe, Below 100°C +/-2°C
Temperature Probe	Available as an Accessory
Calibration	1 Plate/1 Probe Temperature Calibration
Power	220V-240V~ 4A 50/60Hz
