

# MOISTURE ANALYZER – MB120 UK - 30303304

## SPECIFICATIONS

Maximum Capacity	120 g
Readability Moisture Content	1 mg/0.01%
Pan Size	3.5 in (90 mm)
Heater Type	Halogen Lamp
Communication	USB device; USB host x 2; RS232
Dimensions	13.98 in x 7.09 in x 8.31 in (355 mm x 180 mm x 211 mm) (LxHxW)
Display	Full-Color 4.3 in Touchscreen WQVGA
Heating Profiles	Ramp; Standard; Step; Fast

---

In-Use Cover	Included
Moisture Range	0 % ; 100 %
Net Weight	11.53 lb (5.23 kg)
Operating Range (Temp.)	10 °C ; 40 °C
Power	230V
Sample Size	0.5 g ; 120 g
Shut-Off Criteria	Timed; Manual; Custom - User Set Weight Loss/Time; Custom - User Set Moisture% Loss/Time; Automatic - Preset W Loss/Time
Temperature Range	40 °C ; 230 °C
Test Library	100

---

---

Units of Measurement

%DC; Gram; %RG; %MC

---

Working Environment

50°F – 104°F, 85%RH, non-condensing (10°C – 40°C, 85% RH, non-condensing)

---