

# MOISTURE ANALYZER – MB25 UK – 80252472

## SPECIFICATIONS

Maximum Capacity	110 g
Readability Moisture Content	5 mg/0.05%
Pan Size	3.5 in (90 mm)
Heater Type	Halogen Lamp
Communication	RS232 (Included)
Dimensions	11 in x 5 in x 6.5 in (280 mm x 127 mm x 165 mm) (LxHxW)
Display	Backlit LCD
In-Use Cover	Included

---

Net Weight	5 lb (2.1 kg)
Power	AC Power (Included)
Shut-Off Criteria	Timed; Manual; Automatic - Preset Weight Loss/Time
Temperature Range	50 °C – 160 °C
Units of Measurement	Gram
Working Environment	50°F – 104°F, 80%RH, non-condensing (10°C – 40°C, 80%RH, non-condensing)

---