

NAVIGATOR PORTABLE BALANCES – NV3202 – 30615586

SPECIFICATIONS

Maximum Capacity	3,200 g
Readability	0.01 g
Pan Size	5.7 in x 7.5 in (144 mm x 190 mm)
Battery Life	200 Hours with Disposable Batteries
Communication	Ethernet (Available as an Accessory); RS232 (Available as an Accessory); USB Device (Available as an Accessory)
Auxiliary Display Model	Not Applicable
Dimensions	9.1 in x 2.8 in x 8 in (230 mm x 70 mm x 204 mm) (LxHxW)
Display	LCD Display with white backlight

In-Use Cover	No
Legal for Trade	Not applicable
Linearity \pm	0.05 g
Net Weight	2.2 lb (1 kg)
Pan Construction	Stainless Steel
Power	4 C (LR14) Batteries (Not Included); AC Adapter (Included)
Repeatability, typical	0.02 g
Stabilization Time	2 s
Tare Range	To capacity by subtraction
Transportation Case	Available as an Accessory

Units of Measurement

Ounce Troy; Pennyweight; Newton; Gram; Kilogram; Grain; Pound; Ounce; Carat; Pound:Ounce

Working Environment

50°F – 104°F, 85%RH, non-condensing (10°C – 40°C, 85% RH, non-condensing)

