

# PR SERIES PRECISION – PR4202M – 30524611

## SPECIFICATIONS

Maximum Capacity	4,200 g
Readability	0.01 g
Readability (Certified)	0.1 g
Pan Size	7.1 in (180 mm)
Internal Calibration	InCal™ - Semi Automatic
Draftshield	—
Auxiliary Display Model	Available as an Accessory
Battery Life	Not applicable

---

Communication	RS232
Dimensions	12.5 in x 3.7 in x 7.9 in (317 mm x 93 mm x 201 mm) (LxHxW)
Display	LCD with backlight
In-Use Cover	No
Legal for Trade	OIML
Linearity $\pm$	0.02 g
Net Weight	7.7 lb (3.5 kg)
Pan Construction	Stainless Steel
Power	AC Adapter (Included)
Repeatability, typical	0.01 g

---

---

Stabilization Time	1 s
Tare Range	To capacity by subtraction
Units of Measurement	Gram; Kilogram; Carat
Working Environment	50°F – 104°F, 80%RH, non-condensing; 50°F – 86°F, 80%RH, non-condensing (10°C – 40°C, 80%RH, non-condensing; 10°C – 30°C, 80%RH, non-condensing)

---

