

PIONEER PRECISION – PX3202 – 30429807

SPECIFICATION

Maximum Capacity	3,200 g
Readability	0.01 g
Pan Size	180 mm
Internal Calibration	InCal™ - Semi Automatic
Draftshield	—
Auxiliary Display Model	Available as an Accessory
Battery Life	Not applicable
Communication	USB; RS232

Dimensions	321 mm x 98 mm x 209 mm (LxHxW)
Display	2-Line LCD with Backlight
In-Use Cover	Included
Legal for Trade	No
Linearity ±	0.02 g
Minimum Weight (USP, 0.1%, typical)	20 g
Minimum Weight	8.2 g
Net Weight	7.7 lb (3.5 kg)
Pan Construction	Stainless Steel

Power	AC Adapter (Included)
Repeatability, typical	0.01 g
Stabilization Time	1 s
Tare Range	To capacity by subtraction
Units of Measurement	Singapore Tael; Pennyweight; Ounce Troy; Kilogram; Grain; Tical; Custom; Newton; Momme; Taiwan Tael; Baht; Gram; Hong Kong Tael; Pound; Tola; Mesghal; Ounce; Carat
Working Environment	operability guaranteed between 41°F and 104°F; 50°F – 86°F, 80%RH, non-condensing (operability guaranteed between 5°C and 40°C; 10°C – 30°C, 80%RH, non-condensing)

