

# PIONEER SEMI-MICRO – PX125D – 30478196

## SPECIFICATIONS

Maximum Capacity	52 g/120 g
Readability	0.01 mg; 0.1 mg
Readability (Certified)	0.001 g
Pan Size	80 mm
Internal Calibration	AutoCal™ - Automatic
Draftshield	Included
Auxiliary Display Model	Available as an Accessory
Battery Life	Not applicable

---

Communication	USB; RS232
Dimensions	321 mm x 309 mm x 209 mm (LxHxW)
Display	2-Line LCD with Backlight
In-Use Cover	Included
Legal for Trade	OIML
Linearity $\pm$	0.0001 g
Minimum Weight (USP, 0.1%, typical)	20 mg
Net Weight	10 lb (4.5 kg)
Pan Construction	Stainless Steel

---

---

Power	AC Adapter (Included)
Repeatability, typical	0.02 mg; 0.1 mg
Stabilization Time	10 s
Tare Range	To capacity by subtraction
Units of Measurement	Singapore Tael; Ounce Troy; Pennyweight; Grain; Kilogram; Tical; Custom; Milligram; Momme; Newton; Taiwan Tael; Gram; Hong Kong Tael; Pound; Tola; Mesghal; Carat; Ounce
Working Environment	operability guaranteed between 41°F and 104°F; 50°F – 86°F, 80%RH, non-condensing (operability guaranteed between 5°C – 40°C; 10°C – 30°C, 80%RH, non-condensing)

---

