

IMPACT PRINTER – SF40A – 30064202

SPECIFICATIONS

Printer Type	Impact Dot Matrix
Interface Type	RS232
Print Speed	1.2 lines per second (1.2 lines per second)
Print Width	1.89 in (48 mm)
Dimensions	7.99 in x 2.87 in x 4.72 in (203 mm x 73 mm x 120 mm) (LxHxW)
Net Weight	1.55 lb (704 g)
Power	AC Adapter (Included)
Working Environment	32°F – 72°F, 30% – 80% RH, non-condensing (0°C – 40°C, 30% – 80%RH, non-condensing)