

SCOUT PORTABLE BALANCES – STX1202 – 30253063

SPECIFICATIONS

Maximum Capacity	1,200 g
Readability	0.01 g
Pan Size	5.5 in x 6.7 in (140 mm x 170 mm)
Battery Life	6 Hours with Disposable Batteries
Communication	Ethernet (Available as an Accessory); RS232 (Available as an Accessory); USB Host (Available as an Accessory); Bluetooth (Available as an Accessory); USB Device (Available as an Accessory)
Dimensions	8.8 in x 2.1 in x 8 in (224 mm x 54 mm x 202 mm) (LxHxW)
Display	Full-Color 4.3 in Touchscreen WQVGA
In-Use Cover	Available as an Accessory

Legal for Trade	Not Applicable
Linearity \pm	0.03 g
Net Weight	2.2 lb (1 kg)
Pan Construction	Stainless Steel
Power	4 AA (LR6) Batteries (Not Included); AC Adapter (Included)
Repeatability, typical	0.02 g
Stabilization Time	1.5 s
Tare Range	To capacity by subtraction
Transportation Case	Available as an Accessory

Units of Measurement	Ounce Troy; Pennyweight; Newton; Grain; Gram; Kilogram; Pound; Carat; Custom; Ounce; Pound:Ounce
Working Environment	50°F – 104°F, 80%RH, non-condensing (10°C – 40°C, 80%RH, non-condensing)
Auxiliary Display Model	Available as an Accessory

