

SCOUT PORTABLE BALANCES – STX223 – 30253059

SPECIFICATIONS

Maximum Capacity	220 g
Readability	0.001 g
Pan Size	3.7 in (93 mm)
Battery Life	6 Hours with Disposable Batteries
Communication	Ethernet (Available as an Accessory); RS232 (Available as an Accessory); USB Host (Available as an Accessory); Bluetooth (Available as an Accessory); USB Device (Available as an Accessory)
Dimensions	8.7 in x 4.1 in x 8 in (222 mm x 103 mm x 202 mm) (LxHxW)
Display	Full-Color 4.3 in Touchscreen WQVGA
In-Use Cover	Available as an Accessory

Legal for Trade	Not Applicable
Linearity \pm	0.003 g
Net Weight	2.2 lb (1 kg)
Pan Construction	Stainless Steel
Power	4 AA (LR6) Batteries (Not Included); AC Adapter (Included)
Repeatability, typical	0.002 g
Stabilization Time	1.5 s
Tare Range	To capacity by subtraction
Transportation Case	Available as an Accessory

Units of
Measurement

Milligram; Ounce Troy; Pennyweight; Newton; Grain; Gram; Carat; Custom; Ounce

Working
Environment

50°F – 104°F, 80%RH, non-condensing (10°C – 40°C, 80%RH, non-condensing)

Auxiliary Display
Model

Available as an Accessory

