

# VORTEX HEAVY DUTY MIXER – VXHDDG - 30392136

## SPECIFICATIONS

|             |  |
|-------------|--|
| Control     | Digital  |
| Duty Rating | Continuous   |
| Orbit       | 0.2 in (4.9 mm)  |
| Speed Range | On Mode, 300 rpm – 2500 rpm; Touch Mode, 300 rpm – 3500 rpm  |
| Capacity    | 1.1 kg   |
| Timer       | 1 second - 9999 minutes                                      |
| Dimensions  | 9.5 in x 6.3 in x 6.6 in (24.1 cm x 16 cm x 16.8 cm) (LxHxW) |
| Net Weight  | 14.8 lb (6.7 kg)   |

---

|                      |   |
|----------------------|---|
| Power                | 230V, 0.13A, 50/60Hz  |
| Power Consumption    | 30 W  |
| Safety Certification | CE; TUV   |
| Working Environment  | 39°F – 104°F, 20% – 85%RH, non-condensing (4°C – 40°C, 20% – 85%RH, non-condensing) |

---