

VORTEX MULTI TUBE MIXER – VXMTDG - 30392173

SPECIFICATIONS

Control	Digital
Duty Rating	Continuous
Orbit	0.1 in (3.6 mm)
Speed Range	500 rpm – 2500 rpm
Capacity	4.5 kg
Timer	1 second - 9999 minutes
Dimensions	9.5 in x 16 in x 15.1 in (24.1 cm x 40.6 cm x 38.4 cm) (LxHxW)
Net Weight	44.8 lb (20.3 kg)

Power	230V, 1A, 50/60Hz
Power Consumption	100 W
Safety Certification	CE; TUV
Tray Size	7.25 in x 12.25 in (18.4 cm x 31.1 cm) (LxW)
Working Environment	39°F – 104°F, 20% – 85%RH, non-condensing (4°C – 40°C, 20% – 85%RH, non-condensing)
